



Disposable & Reusable Electrosurgical Pencil




- ☐ Standard features that compatible with most \varnothing 2.36mm electrodes and electrosurgical units.
- ☐ Flexible 3m or 5m cord for easy set-up and drapeability.
- ☐ Ergonomic design for comfortable grip and operation in surgery.
- ☐ All monopolar ES pencil rated accessory voltage of 5kVp.
- ☐ Disposable types in sterile packs and Reusable types are available.

ES Pencil




Disposable

- PD702T (BOX/10PCS)
Push Button Pencil with Teflon coated blade
Cable length: 3 meter
Plug type: 3 prong plug
Fit for 2.36mm Electrode
Sterile packing
- PD702B (BOX/10PCS)
Push Button Pencil with Standard blade




Disposable

- PD700T (BOX/10PCS)
Push Button Pencil with Teflon coated blade
Cable length: 3 meter
Plug type: 3 prong plug
Fit for 2.36mm Electrode
Sterile packing
- PD700B (BOX/10PCS)
Push Button Pencil with Standard blade



Disposable

- **Push Button Pencil with Stamping blade**
Plug type: 3 prong plug
Fit for 2.36mm Electrode
Sterile packing
- PD513 (BOX/10PCS)
Cable length: 3 meter
- PD513P5 (BOX/10PCS)
Cable length: 5 meter




Disposable

- PD513T (BOX/10PCS)
Push Button Pencil with Teflon coated blade
Cable length: 3 meter
Plug type: 3 prong plug
Fit for 2.36mm Electrode
Sterile packing




Reusable

- PR306P (Box/10pcs)
Push Button Pencil with Reusable blade
Reuse 30 times
Cable length: 3 meter
Plug type: 3 prong plug
Sterile packing



- **3 prong plug included with ES pencil**
Cable length: 3 or 5 meter



- **Tip Polisher**
• TC50501
Size: 5 × 5 cm