


**Upspring**



# **PRODUCT CATALOG**

**NEUROSURGERY SOLUTION**

PRECISION | INNOVATION | TRUST

# Upspring Medical

At Upspring Medical, we are committed to delivering cutting-edge surgical solutions to the Turkish & global healthcare sector. Established in Ankara in 2024, our office is officially registered as a Medical Sales Center under the Sincan District Health Directorate and listed in the ÜTS system. Following full regulatory approval in April 2025, we launched nationwide operations, combining global manufacturing excellence with local expertise to provide reliable, high-performance solutions.

Our products are manufactured through a strategic partnership with Promed Supply Co., Ltd., a Hong Kong-based leader in electrosurgical innovation since 1997. Promed operates ISO 13485 and CE-certified facilities in China, ensuring top-tier quality, full compliance, and global production standards.

**Precision.  
Innovation.  
Trust**

*Delivering Globally Reliable  
Electrosurgical Solutions*

## Local Commitment with Global Expertise

At Upspring Medical, we leverage our global manufacturing expertise to deliver exceptional quality and service to the Turkish & global market. By combining Promed's international production strength with localized support, we offer tailored solutions and responsive service to meet the needs of Turkish & global distributors and healthcare professionals.

## End-to-End Manufacturing & Reliable Supply Chain

Promed supply oversees the entire production process from raw materials to final sterilized packaging ensuring **full traceability**, **consistent quality**, and **on-time delivery**.





## Comprehensive Electrosurgical Solutions

- ◆ Electrosurgical Pencils (Disposable & Reusable) & Bipolar Forceps
- ◆ Specialty Electrodes (General, LLETZ, Micro, Arthroscopic) & Grounding Pads & Cables
- ◆ Laparoscopic Instruments & Electrodes & Skin Staplers & Removers
- ◆ Smoke Evacuation Systems (Pencils and Forceps)

These solutions support a wide range of surgical specialties

- ◆ From **General Surgery** to **Gynecology, Aesthetic**, and **Abdominal procedures**
- ◆ Ensuring precision, safety, and performance across every operation.

# TABLE OF CONTENTS



## NEUROSURGERY SOLUTION

*High-precision and low-damage tools for neurological operations*

- 01** Disposable & Reusable Electrosurgical Pencil
- 05** Reusable Bipolar Forceps
- 09** Disposable Suction Coagulators
- 10** Disposable Suction Irrigator
- 11** Disposable Ground Pad & Cable
- 15** Disposable Skin Stapler
- 16** Disposable Skin Stapler Remover
- 17** Disposable Smoke Evacuation Pencil
- 19** Disposable Smoke Evacuation Forceps





# Disposable Electrosurgical Pencil

*Power At Your Fingertips*



- Equipped with stainless steel and Teflon-coated blade to reduce tissue adhesion.
- Lightweight, ergonomic design with push-button activation for precise control.
- Universally compatible with  $\varnothing 2.36\text{mm}$  electrodes and most ESU generators.
- Offered with 3m or 5m flexible cable options for enhanced maneuverability.
- CE-certified and rated for up to 5 kVp accessory voltage.



## Disposable Electrosurgical Pencil

- Fit for with  $\varnothing$  2.36mm Electrode (Stainless Steel / Teflon Coated Blade).
- Sterile Packing
- Packing:10pcs/Box



PRODUCT CODE	BLADE TYPE	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
<b>PD702T</b>	Teflon Coated Blade	168MM	3M	3 Prong Plug
<b>PD702B</b>	Standard Blade	168MM	3M	3 Prong Plug



PRODUCT CODE	BLADE TYPE	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
<b>PD400T</b>	Teflon Coated Blade	168MM	3M	3 Prong Plug
<b>PD400S</b>	Standard Blade	168MM	3M	3 Prong Plug



PRODUCT CODE	BLADE TYPE	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
<b>PD513</b>	Stamping Blade	168MM	3M	3 Prong Plug
<b>PD513P5</b>	Stamping Blade	168MM	5M	3 Prong Plug



PRODUCT CODE	BLADE TYPE	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
<b>PD513T</b>	Teflon Coated Blade	168MM	3M	3 Prong Plug

# Disposable Microfin Electrosurgical Pencil

Exceptional Precision for Micro-Electrosurgery



Cut and coagulation buttons are conveniently positioned near the surgical site for improved control, precision, and smooth electrode handling.

Compatible with ultra-fine  $\varnothing 1.6\text{mm}$  electrode for superior surgical precision.

Lightweight, compact design for easy handling and reduced fatigue.

**PD833W -**  
 **$\varnothing 1.6\text{MM}$**

**PD833 -**  
 **$\varnothing 2.36\text{MM}$**

Compatible with standard  $\varnothing 2.36\text{mm}$  electrodes for versatile surgical applications.

\*\*\* Electrode is not included with the ESU pencil \*\*\*

## Disposable Microfin Electrosurgical Pencil

- Sterile Packing
- Packing: 10pcs/Box



PRODUCT CODE	ELECTRODE DIAMETER	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
PD833W	$\varnothing 1.6\text{MM}$	106MM	3M	3 Prong Plug
PD833	$\varnothing 2.36\text{MM}$	106MM	3M	3 Prong Plug



# Reusable Electrosurgical Pencil

Reliable, Reusable, and Cost-Effective



NEUROSURGERY SOLUTION

Reusable for up to 30/50/100 cycles, reducing long-term operating costs.

Push-button activation enables precise and efficient electrosurgical performance.


Ergonomic and lightweight design for enhanced surgeon comfort and control.

Compatible with 2.36mm electrodes and most standard electrosurgical units.

## Reusable Electrosurgical Pencil

- Fit with  $\varnothing$  2.36mm Electrode
- Sterile Packing
- Packing: 10pcs/Box

134°C  
Autoclave



PRODUCT CODE	BLADE TYPE	REUSABLE	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
PR306P	Standard Reusable Blade	30 Cycles	163MM	3M	3 Prong Plug
PR307P	Standard Reusable Blade	50 Cycles	163MM	3M	3 Prong Plug
PR310P	Standard Reusable Blade	100 Cycles	163MM	3M	3 Prong Plug



# Reusable Bipolar Forceps

*Durable Precision For Bipolar Coagulation*



**Promed's Reusable Bipolar Forceps** are available in a variety of lengths, tip sizes, and designs, including Bayonet, Jewelers, Adson, and Straight, to suit diverse surgical needs.

Engineered for optimal performance, the polished stainless-steel tips reduce eschar buildup, ensuring smooth operation. The ergonomic, balanced handles minimize hand fatigue, enhancing comfort and precision. Designed for durability, they withstand up to 20 autoclave cycles.



## Polished Stainless Steel Tips

Secure grasping and coagulation at lower power, reducing tissue damage and eschar buildup.



## Insulated Stainless Steel Body

Prevents electrical leakage for enhanced safety.



## Precise Handle Spring Tension

Ensuring consistent dissection pressure and maintaining accurate tip alignment.




## Platform-Like Tip Design

Increases resistance and enhances coagulation with controlled heat.



## Ergonomic, Well-Balanced Handles


Provides balance and comfort reducing hand fatigue.





STYLE	LENGTH(MM)	TIP WIDTH (MM)
Bayonet Straight	180/200/220/240MM	0.5/1.0/2.0MM
Bayonet Angled Down	200/220/240MM	0.5/1.0/2.0MM
Bayonet Angled Up	200/220/240MM	0.5/1.0/2.0MM
Straight	140/180/200/220/240MM	0.5/1.0/2.0MM
Adson	130MM	0.5/1.0MM
Jeweler	110MM	0.5MM

# Straight Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non-sterile, 1 each per pack





FST0205, Length 200mm, Straight, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
<b>Straight Tip, 2 Pins Connector</b>					
<b>0.5MM</b>	FST0145	FST0185	FST0205	/	/
<b>1.0MM</b>	FST0141	FST0181	FST0201	FST0221	FST0241
<b>2.0MM</b>	FST0142	FST0182	FST0202	FST0222	FST0242
<b>Straight Tip, Flat Connector</b>					
<b>0.5MM</b>	FST1145	FST1185	FST1205	/	/
<b>1.0MM</b>	FST1141	FST1181	FST1201	FST1221	FST1241
<b>2.0MM</b>	FST1142	FST1182	FST1202	FST1222	FST1242



FAG0205, Length 200mm, Straight, Angled Tip, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
<b>Straight, Angled Tip, 2 Pins Connector</b>					
<b>0.5MM</b>	/	FAG0185	FAG0205	/	/
<b>1.0MM</b>	/	FAG0181	FAG0201	FAG0221	FAG0241
<b>2.0MM</b>	/	FAG0182	FAG0202	FAG0222	FAG0242
<b>Straight, Angled Tip, Flat Connector</b>					
<b>0.5MM</b>	/	FAG1185	FAG1205	/	/
<b>1.0MM</b>	/	FAG1181	FAG1201	FAG1221	FAG1241
<b>2.0MM</b>	/	FAG1182	FAG1202	FAG1222	FAG1242

# Bayonet Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack



FBD0205, Length 200mm, Bayonet, Angled Down, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
<b>Bayonet, Angled Down, 2 Pins Connector</b>					
<b>0.5MM</b>	/	/	FBD0205	FBD0225	/
<b>1.0MM</b>	/	/	FBD0201	FBD0221	FBD0241
<b>2.0MM</b>	/	/	FBD0202	FBD0222	/
<b>Bayonet, Angled Down, Flat Connector</b>					
<b>0.5MM</b>	/	/	FBD1205	FBD1225	/
<b>1.0MM</b>	/	/	FBD1201	FBD1221	FBD1241
<b>2.0MM</b>	/	/	FBD1202	FBD1222	/





FBS0205, Length 200mm, Bayonet, Straight, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
<b>Bayonet, Straight Tip, 2 Pins Connector</b>					
<b>0.5MM</b>	/	FBS0185	FBS0205	FBS0225	FBS0245
<b>1.0MM</b>	/	FBS0181	FBS0201	FBS0221	FBS0241
<b>2.0MM</b>	/	FBS0182	FBS0202	FBS0222	/
<b>Bayonet, Straight Tip, Flat Connector</b>					
<b>0.5MM</b>	/	FBS1185	FBS1205	FBS1225	FBS1245
<b>1.0MM</b>	/	FBS1181	FBS1201	FBS1221	FBS1241
<b>2.0MM</b>	/	FBS1182	FBS1202	FBS1222	/




FBU0205, Length 200mm, Bayonet, Angled Up, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
<b>Bayonet, Angled Up, 2 Pins Connector</b>					
<b>0.5MM</b>	/	/	FBU0205	FBU0225	/
<b>1.0MM</b>	/	/	FBU0201	FBU0221	FBU0241
<b>2.0MM</b>	/	/	FBU0202	FBU0222	/
<b>Bayonet, Angled Up, Flat Connector</b>					
<b>0.5MM</b>	/	/	FBU1205	FBU1225	/
<b>1.0MM</b>	/	/	FBU1201	FBU1221	FBU1241
<b>2.0MM</b>	/	/	FBU1202	FBU1222	/

## Jewelers

### Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack




MJE0115, Length 110mm, Jewelers, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)
	110MM
<b>Jeweler, 2 Pins Connector</b>	
<b>0.5MM</b>	MJE0115
<b>1.0MM</b>	/
<b>Jeweler, Flat Connector</b>	
<b>0.5MM</b>	MJE1115
<b>1.0MM</b>	/

## Adson

### Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack



MAD0135, Length 130mm, Adson, 0.5mm Tip Width, 2 Pins Connector


TIP WIDTH (MM)	FORCEPS LENGTH(MM)
	130MM
<b>Adson, 2 Pins Connector</b>	
<b>0.5MM</b>	MAD0135
<b>1.0MM</b>	MAD0131
<b>Adson, Flat Connector</b>	
<b>0.5MM</b>	MAD1135
<b>1.0MM</b>	MAD1131

# Suction Coagulator

Dual Function :Suction & Coagulation




## Tube Diameter 8Fr/10Fr/12Fr




<b>Type</b>	Footswitch or Hand Switch
<b>Size Available</b>	8Fr,10Fr,12Fr
<b>Suction:</b>	Place Finger Over Suction Vent on Handle
<b>Handle:</b>	Ergonomically Designed White ABS
<b>Cannula Insulation:</b>	Blue Polyolefin
<b>Obturator:</b>	Polypropylene Cleaning Tube and Protective Insert for Bending
<b>Suction Tubing:</b>	6mm or 5mm Tubing Compatible
<b>Cable:</b>	3m(10ft)/5m(15ft)for both Footswitch and Handswitch

## Handswitch Suction Coagulators

PRODUCT	CODE	TUBE DIAMETER	CABLE LENGTH(M)	PACKING
	<b>SC08H</b>	8fr	3M	15/Case
	<b>SC10H</b>	10fr	3M	15/Case
	<b>SC12H</b>	12fr	3M	15/Case

## Footswitch Suction Coagulators

PRODUCT	CODE	TUBE DIAMETER	CABLE LENGTH(M)	PACKING
	<b>SC08F</b>	8fr	3M	15/Case
	<b>SC10F</b>	10fr	3M	15/Case
	<b>SC12F</b>	12fr	3M	15/Case






# Disposable Suction Irrigator

Clean Precision , Seamless Control

Engineered to streamline surgical workflows, Promed's Suction Irrigator integrates simultaneous suction and irrigation in a unified system. Featuring an ergonomic design and standardized connectivity, it maintains optimal visibility during procedures while supporting precision in fluid evacuation and targeted irrigation.



NEUROSURGERY SOLUTION



## Order Information


PRODUCT CODE	PRODUCT DESCRIPTION	PACKING
SI0528PD	Disposable Suction Irrigator, $\varnothing$ 5.36mmx28cm Composites Probe, Double Spikes	20/Case
SI0533PD	Disposable Suction Irrigator, $\varnothing$ 5.36mmx 33cm Composites Probe, Double Spikes	20/Case
S10545PD	Disposable Suction Irrigator, $\varnothing$ 5.36mmx45cm Composites Probe, Double Spikes	20/Case
S10528PS	Disposable Suction Irrigator, $\varnothing$ 5.36mm X28cm Composites Probe, Single Spike	20/Case
S10533PS	Disposable Suction Irrigator, $\varnothing$ 5.36mmx 33cm Composites Probe, Single Spike	20/Case
SI0545PS	Disposable Suction Irrigator, $\varnothing$ 5.36mm X45cm Composites Probe, Single Spike	20/Case

# Disposable Ground Pad & Cable

*Skin-Friendly · Single-Use · Burn-Reducing*



Promed's Disposable Grounding Pad offers safe, reliable grounding during electrosurgery. Its soft foam and hydrogel ensure secure, comfortable adhesion while optimizing energy transfer. Skin-friendly, single-use, and burn risk-reducing.




The **Round edge design** evenly distributes current, reducing concentrated areas and minimizing the risk of patient

The **Biocompatible hydrogel** adhesive provides instant, consistent and comfortable impedance upon skin contact.

The **highly flexible adhesive foam** border and conductive gel ensure secure electrode placement.

**Highly compatible** with most electrosurgical generators, supporting bipolar and monopolar designs in various sizes.



# Disposable Grounding Pad Without Cable

- Non-Sterile Packing
- Packing:1PC/1 Aluminum Pouch, 25pcs/Box




PRODUCT	PRODUCT CODE	STYLE	SHAPE	USE FOR	SIZE (LENGTH X WIDTH) (MM)
	GM100	Monopolar	Vertical	Pediatric	148MM*90MM
	GM900	Monopolar	Vertical	Adult	204MM*106MM
	GB100	Bipolar	Vertical	Pediatric	148MM*90MM
	GB900	Bipolar	Vertical	Adult	204MM*106MM
	GM100H	Monopolar	Horizontal	Pediatric	166MM*90MM
	GM900H	Monopolar	Horizontal	Adult	206MM*104MM
	GB100H	Bipolar	Horizontal	Pediatric	166MM*90MM
	GB900H	Bipolar	Horizontal	Adult	206MM*104MM




# Disposable Grounding Pad Integrated With Cable

- Prewired with single-use cable permanently attached
- Cable Length: 2.8 Meter
- Non-Sterile Packing
- Packing: 1PC/1 Aluminum Pouch, 25pcs/Box



GB910, Disposable Ground Pad, Bipolar, Vertical, Adult

PRODUCT	CODE	STYLE	SHAPE	USE FOR	SIZE (LENGTH X WIDTH) (MM)	PLUG TYPE
	GM110	Monopolar	Vertical	Pediatric	148MM*90MM	Non-REM Plug
	GM910	Monopolar	Vertical	Adult	204MM*106MM	Non-REM Plug
	GB110	Bipolar	Vertical	Pediatric	148MM*90MM	REM- Plug
	GB910	Bipolar	Vertical	Adult	204MM*106MM	REM- Plug
	GM110H	Monopolar	Horizontal	Pediatric	166MM*90MM	Non-REM Plug
	GM910H	Monopolar	Horizontal	Adult	206MM*104MM	Non-REM Plug
	GB110H	Bipolar	Horizontal	Pediatric	166MM*90MM	REM- Plug
	GB910H	Bipolar	Horizontal	Adult	206MM*104MM	REM- Plug



## Disposable Detachable Grounding Pad Cable

- Non-Sterile Packing
- PE Bag, 50pcs/Box



Non-Rem Plug






Rem Plug



Hi-Fi Plug



Clip size  
L×W×H=54×31×11 MM

PRODUCT	CODE	PLUG TYPE	CABLE TYPE	CABLE LENGTH(M)	CLIP SIZE L×W×H (MM)
	<b>GC2VR</b>	Valley Lab REM Plug	PVC	3 Meters	54MM×31MM×11MM
	<b>GC2VN</b>	Valley Lab Non-REM Plug	PVC	3 Meters	54MM×31MM×11MM
	<b>GC2HF6</b>	6.3mm Hi-Fi Plug	PVC	3 Meters	54MM×31MM×11MM



# Disposable Skin Stapler & Remover

Fast & Safe & Consistent




## Disposable Skin Stapler

- Precise staple forming ensures optimal wound approximation and healing.
- Consistent, smooth, and safe actuation provides accurate staple depth.
- Preloaded with medical-grade stainless steel staples for reliable closure.
- EO is sterilized and designed for single-use to guarantee maximum hygiene.
- A transparent window allows easy monitoring of remaining staples.



### Staple Dimensions

<b>Open</b>	<b>14mm</b>  3.2mm
<b>Closed</b>	<b>W 6.5+0.5/-0.2mm</b>  <b>H 4.7+0.5/-0.4mm</b>

PRODUCT CODE	PRODUCT DESCRIPTION	STAPLE	STERILE PACKING
<b>SW904-15W</b>	Disposable Skin Stapler SW904,15W	15PCS	6/Case
<b>SW904-25W</b>	Disposable Skin Stapler SW904,25W	25PCS	6/Case
<b>SW904-35W</b>	Disposable Skin Stapler SW904,35W	35PCS	6/Case

# Disposable Skin Staple Remover

Promed's Disposable Skin Stapler Remover features an angled tip for easy access and precise staple removal. Designed for gentle tissue handling, it minimizes patient discomfort during procedures. Its simple, user-friendly mechanism ensures efficient and safe use by medical professionals. Sterile and disposable, it supports optimal hygiene and convenience in clinical settings.



Precision-Engineered Tip	Designed to securely grip and remove staples without causing trauma to surrounding tissue.
Angled Tip Design	Allows for easy access and optimal visibility, especially in hard-to-reach areas.
Ergonomic Handle	Comfortable, non-slip grip for controlled and accurate removal.
Sterile & Single-Use	Pre-sterilized and disposable to prevent cross-contamination and ensure patient safety.
Lightweight & Easy to Use	Compact and lightweight design minimizes hand fatigue and enhances ease of use during procedures.
Durable Medical-Grade Material	Made from high-quality plastic and stainless steel components for reliable performance.

PRODUCT CODE	PRODUCT DESCRIPTION	STAPLE	STERILE PACKING
SR00D1	Disposable Skin Staple Remover	15PCS	10/Case

# Telescopic Smoke Evacuation Pencil

Slide-To-Extend Technology For Complete Surgical Control



## SmokeFin

Smart · Clear · Ergonomic · Disposable

- **Enhanced Ergonomics**  
Lightweight with 360° swivel and flexible tubing.

- **Slim, Oval Triangular Grip**  
Ensures improved access and visibility.

- **Rocker Switch with Sliding Button**  
Telescopic design that slides along the housing for quick adjustments.


- **Efficiency Smoke Clearance**  
Rapid removal of surgical smoke to maintain a clear operating field.

- **Slide-to-Extend up to 60mm**  
By pushing the button forward, the housing can extend up to 60mm to reach deeper into surgical cavities.


- **Adjustable Length**  
The distance between the sliding button and electrode tip smoothly extends from 68MM to 198mm.

68MM-198MM

Hassle-free  
work-length adjust



Full extension



Pushing the button forward extends the electrode tip, while pulling it backward shortens the reach for precise control.

## Adjustable Reach, Efficient Smoke Evacuation

The **SmokeFin Telescopic Smoke Evacuation Pencil** features a sliding button that enables surgeons to instantly adjust the working length from 68mm to 198mm without changing the electrode. Its innovative telescopic mechanism ensures seamless adaptation during procedures, providing superior visibility, enhanced ergonomics, and effortless control.



With common ES Pencil

With SmokeFin

## Accessories

### Nose Cone



- 3 sizes of nose cones fit for 63.5/69/120MM electrodes
- Transparent, equipped with radiopaque identifier (xray-detectable)

### Electrodes



- Standard 50MM  
Optional 63.5/69/120MM
- Modified blade electrode (shrink wrap insulation), coated or non-coated

### Cable and Tubing



- Standard 3 meter or 4.5 meter
- 22MM filter port adaptor
- Y-connection


### LLETZ/ LEEP electrode



- Compatible with a range of electrode types to suit different surgical procedures

## Smoke Evacuation Pencil

- Sterile Packing
- Packing:5pcs/Box



PRODUCT CODE	HOUSING LENGTH	ELECTRODE LENGTH	BLADE TYPE	CABLE LENGTH(M)	PLUG TYPE
<b>PS300</b>	195MM	50MM	Non-Coated	3M	3 Prong Plug
<b>PS300T</b>	195MM	50MM	Non-Sticked	3M	3 Prong Plug
<b>PS400</b>	195MM	63.5MM	Non-Coated	3M	3 Prong Plug
<b>PS400T</b>	195MM	63.5MM	Non-Sticked	3M	3 Prong Plug
<b>PS500</b>	195MM	69MM	Non-Coated	3M	3 Prong Plug
<b>PS500T</b>	195MM	69MM	Non-Sticked	3M	3 Prong Plug
<b>PS600</b>	195MM	120MM	Non-Coated	3M	3 Prong Plug
<b>PS600T</b>	195MM	120MM	Non-Sticked	3M	3 Prong Plug



# Disposable Smoke Evacuation Forceps

Precision Grasping With Smoke Evacuation



**2 - In - 1 Solution**  
for electrosurgery & smoke evacuation  
reduces surgical plume

Hand Switch

360° Swivel

*Monopolar*

*Bipolar*

Smoke Plume Capture

Available In Monopolar  
Non-Evacuation version  
with Non-coated Tip



## Accessories

### ▪ Tips Size



1.20mm  
Serrated




1.20mm  
Smooth



0.75mm  
Serrated



0.75mm  
Smooth




0.5mm  
Smooth



0.1mm  
Smooth

### ▪ Compatible Accessories

#### Cable and Tubing



- Cable and tubing available in 3 meters and 4.5 meters
- Tubing and connector available in different sizes

#### Plug



Monopolar




Bipolar

- Compatible with all the generator such as ValleyLab, Bovie, and Conmed



## Features



### Available in Both Bipolar and Monopolar Configurations

Both Bipolar and Monopolar configurations are available catering to different needs.



### Smoke Evacuation

Extracting hazardous surgical smoke at the tip of the forceps with incorporated smoke tube.

*Non-Evacuation monopolar forceps is also available.*



### Hand Activation Control


Completely hand-activated control by squeezing and holding the buttons on the sides;

Footswitch is not necessary during surgery.



### Flexible Tubing and Lightweight

Reducing the operator fatigue while maximizing the flexibility by the integration of smoke tubing with 360° swivel connection.



### Tips with Optional PTFE Coating

Reducing the risk of building up eschar on the tips and delivering reliable electrosurgical performance.

*PTFE coating is available in Monopolar Forceps to eliminate tissue adhesion further.*

## Disposable Smoke Evacuation Forceps



PRODUCT CODE	PRODUCT DESCRIPTION	TIP TYPE	TIP SIZE	CABLE LENGTH (M)	PLUG TYPE
<b>FBS001</b>	Single Use Hand Switch Bipolar Forceps with Smoke Evacuation	Non-coated	0.1MM Smooth Tip	3M	Bipolar
<b>FBS005</b>	Single Use Hand Switch Bipolar Forceps with Smoke Evacuation	Non-coated	0.5MM Smooth Tip	3M	Bipolar
<b>FBS006</b>	Single Use Hand Switch Bipolar Forceps with Smoke Evacuation	Non-coated	0.75MM Smooth Tip	3M	Bipolar
<b>FBS007</b>	Single Use Hand Switch Bipolar Forceps with Smoke Evacuation	Non-coated	1.2MM Smooth Tip	3M	Bipolar
<b>FMS001</b>	Single Use Hand Switch Monopolar Forceps with Smoke Evacuation	Teflon coated	0.1MM Smooth Tip	3M	Monopolar
<b>FMS005</b>	Single Use Hand Switch Monopolar Forceps with Smoke Evacuation	Teflon coated	0.5MM Smooth Tip	3M	Monopolar
<b>FMS006</b>	Single Use Hand Switch Monopolar Forceps with Smoke Evacuation	Teflon coated	0.75MM Smooth Tip	3M	Monopolar
<b>FMS007</b>	Single Use Hand Switch Monopolar Forceps with Smoke Evacuation	Teflon coated	1.2MM Smooth Tip	3M	Monopolar



## **Upspring BRI Corporation Medical**

Tandoğan Mahallesi Billur  
caddesi Özer İş Hanı No:  
12/15 06930 Sincan Ankara, Turkey

Email:

[sales@upspringpurification.com](mailto:sales@upspringpurification.com)

Website:

[www.upspringmedical.com](http://www.upspringmedical.com)