

# Disposable LLETZ Electrode

Optimized For Cervical Procedure



Promed's gynecology loop electrode ensures precise excision with minimal thermal damage, promoting faster healing and reducing scarring. It delivers consistent energy for accurate cuts, minimizing tissue adhesion and improving procedural efficiency.

## DISPOSABLE LLETZ ELECTRODE

|                               |                                       |
|-------------------------------|---------------------------------------|
| Shaft Diameter                | 2.36MM                                |
| Material                      | Tungsten Wire + Medical-Grade Polymer |
| Design                        | Overmold Hex Lock                     |
| Shape Avalale                 | Loop & Square Shape                   |
| Loop Electrode Diameter(MM)   | 5MM,10MM,15MM,20MM                    |
| Square Electrode Diameter(MM) | 5MM,8MM,10MM                          |
| Packing                       | 1pcs/Pouch,10pcs/Box                  |



## FEATURES

|                                        |                                                                                                                        |
|----------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Fine Tungsten Wire Construction</b> | High conductivity and durability for precise cutting, minimizing adhesion and ensuring consistent energy.              |
| <b>Precision Design</b>                | 2.36mm shaft for optimal maneuverability; sterile and single-use to reduce contamination risks.                        |
| <b>Overmold Shaft with HEX Lock</b>    | Ergonomic, non-slip grip for enhanced control during procedures                                                        |
| <b>Multiple Sizes &amp; Shapes</b>     | Available in various loop diameters (e.g., 5mm, 10mm, 15mm, 20mm) for excision, and square electrodes for coagulation. |



EL2V64, Loop (A)5MM x (B)5MM

| PRODUCT CODE | PRODUCT DESCRIPTION         | SHAFT DIAMETER |
|--------------|-----------------------------|----------------|
| EL2V64       | LLETZ Loop 5*5 MM Electrode | Ø2.36MM        |



EL2V60, Loop (A)10MM x (B)7MM

| PRODUCT CODE | PRODUCT DESCRIPTION          | SHAFT DIAMETER |
|--------------|------------------------------|----------------|
| EL2V61       | LLETZ Loop 10*5MM Electrode  | Ø2.36MM        |
| EL2V60       | LLETZ Loop 10*7MM Electrode  | Ø2.36MM        |
| EL2V65       | LLETZ Loop 10*10MM Electrode | Ø2.36MM        |



EL2V6B, Loop (A)15MM x (B)5MM

| PRODUCT CODE | PRODUCT DESCRIPTION          | SHAFT DIAMETER |
|--------------|------------------------------|----------------|
| EL2V6B       | LLETZ Loop 15*5MM Electrode  | Ø2.36MM        |
| EL2V6C       | LLETZ Loop 15*8MM Electrode  | Ø2.36MM        |
| EL2V6P       | LLETZ Loop 15*10MM Electrode | Ø2.36MM        |
| EL2V6J       | LLETZ Loop 15*12MM Electrode | Ø2.36MM        |
| EL2V6K       | LLETZ Loop 15*15MM Electrode | Ø2.36MM        |



EL2V6F, Loop (A)20MM x (B)15MM

| PRODUCT CODE | PRODUCT DESCRIPTION          | SHAFT DIAMETER |
|--------------|------------------------------|----------------|
| EL2V6D       | LLETZ Loop 20*8MM Electrode  | Ø2.36MM        |
| EL2V6E       | LLETZ Loop 20*10MM Electrode | Ø2.36MM        |
| EL2V6H       | LLETZ Loop 20*12MM Electrode | Ø2.36MM        |
| EL2V6F       | LLETZ Loop 20*15MM Electrode | Ø2.36MM        |
| EL2V6G       | LLETZ Loop 20*20MM Electrode | Ø2.36MM        |



EL2V74, Square (A)5MM x (B)5MM

| PRODUCT CODE | PRODUCT DESCRIPTION          | SHAFT DIAMETER |
|--------------|------------------------------|----------------|
| EL2V74       | LLETZ Square 5*5MM Electrode | Ø2.36MM        |



EL2V7A, Loop Square (A)10MM x (B)8MM

| PRODUCT CODE | PRODUCT DESCRIPTION            | SHAFT DIAMETER |
|--------------|--------------------------------|----------------|
| EL2V71       | LLETZ Square 10*5MM Electrode  | Ø2.36MM        |
| EL2V7A       | LLETZ Square 10*8MM Electrode  | Ø2.36MM        |
| EL2V75       | LLETZ Square 10*10MM Electrode | Ø2.36MM        |

