

General Disposable Electrode

Safe & Efficient Energy Transmission



Promed's disposable electrosurgical electrodes—offered in blade, needle, and ball shapes—deliver precise cutting and coagulation. Available in various lengths with stainless steel or Teflon-coated shafts, and in over-molded or shrink-wrapped designs, they are fully compatible with standard monopolar pencils, ensuring accuracy, control, and dependable performance in every



Overmold with Hex without Hex collar




Hex Collar with Shrink Wrap



Hex Collar with Shrink Wrap (High Insulation)



Shrink Wrap



Shrink Wrap (High Insulation)



Optional Finishing for Disposable Electrode



Stainless Steel



Teflon Coated

Disposable Blade Electrode

- $\varnothing$ 2.36mm Shaft Diameter
- Packaged Sterile, 50pcs per Box



		PRODUCT CODE	
PRODUCT DESCRIPTION	LENGTH (MM)	FINISHING-STAINLESS STEEL	OPTIONAL FINISHING-TEFLON COATING
Overmold with Hex, Standard Blade			
	69MM	EB22M2	EB22M2T
	153MM	EB26M2	EB26M2T
Hex Collar with Shrink Wrap			
	102MM	EB24P2	EB24P2T
	152MM	EB26P2	EB26P2T
Hex Collar with Shrink Wrap, High Insulation			
	69MM	EB22S2	EB22S2T
	102MM	EB24S2	EB24S2T
	152MM	EB26S2	EB26S2T
Shrink Wrap			
	69MM	EB22M2-H	EB22M2T-H
	102MM	EB24M2-H	EB24M2T-H
	152MM	EB26M2-H	EB26M2T-H
Shrink Wrap, High Insulation			
	69MM	EB22S2-H	EB22S2T-H
	102MM	EB24S2-H	EB24S2T-H
	152MM	EB26S2-H	EB26S2T-H

Disposable Needle Electrode

- ∅ 2.36mm Shaft Diameter
- Packaged Sterile, 50pcs per Box



		PRODUCT CODE	
PRODUCT DESCRIPTION	LENGTH (MM)	FINISHING-STAINLESS STEEL	OPTIONAL FINISHING-TEFLON COATING
Overmold with Hex, Standard Needle			
	69MM	EN22M2	EN22M2T
	153MM	EN26M2	EN26M2T
Hex Collar with Shrink Wrap			
	102MM	EN24P2	EN24P2T
	152MM	EN26P2	EN26P2T
Hex Collar with Shrink Wrap, High Insulation			
	69MM	EN22S2	EN22S2T
	102MM	EN24S2	EN24S2T
	152MM	EN26S2	EN26S2T
Shrink Wrap			
	69MM	EN22M2-H	EN22M2T-H
	102MM	EN24M2-H	EN24M2T-H
	152MM	EN26M2-H	EN26M2T-H
Shrink Wrap, High Insulation			
	69MM	EN22S2-H	EN22S2T-H
	102MM	EN24S2-H	EN24S2T-H
	152MM	EN26S2-H	EN26S2T-H

Disposable Ball Electrode

- $\varnothing$ 2.36mm Shaft Diameter
- Packaged Sterile, 50pcs per Box



			PRODUCT CODE	
PRODUCT DESCRIPTION	LENGTH (MM)	LENGTH (MM)	FINISHING-STAINLESS STEEL	OPTIONAL FINISHING-TEFLON COATING
Overmold with Hex, Standard Ball				
	3MM	69MM	ES22M3	ES22M3T
	5MM	71MM	ES22M4	ES22M4T
	3MM	133MM	ES26M3	ES26M3T
	5MM	136MM	ES26M4	ES26M4T
Hex Collar with Shrink Wrap				
	3MM	69MM	ES22S3	ES22S3T
	5MM	71MM	ES22S4	ES22S4T
	3MM	102MM	ES24S3	ES24S3T
	5MM	107MM	ES24S4	ES24S4T
	3MM	134MM	ES26S3	ES26S3T
	5MM	136MM	ES26S4	ES26S4T
Shrink Wrap				
	3MM	69MM	ES22M3-H	ES22M3T-H
	5MM	71MM	ES22M4-H	ES22M4T-H
	3MM	101MM	ES24M3-H	ES24M3T-H
	5MM	107MM	ES24M4-H	ES24M4T-H
	3MM	134MM	ES26M3-H	ES26M3T-H
	5MM	136MM	ES26M4-H	ES26M4T-H

Microfin Disposable Electrode

Micro-Level Accuracy For Delicate Procedures



The Disposable Microfin Electrode is engineered for precision and performance across multiple specialties, including Aesthetic & Dermatologic Surgery (Skin Lesion Removal, Scar Revision), ENT (Nasal Polyp Excision), and Dentistry (Gingival Sculpting) etc.

Its ultra-fine design allows for exceptional control in delicate procedures, while broad compatibility with electrosurgical systems ensures minimal tissue damage and reduced thermal spread - enhancing both safety and surgical outcomes.

MicroFin



	Round Loop Electrode Available Size • Ø 4MM • Ø 7MM • Ø 10MM • Ø 13MM		Triangle Loop Electrode Available Size • Ø 4MM • Ø 7MM • Ø 10MM		Diamond Loop Electrode Available Size • Ø 3MM • Ø 5MM • Ø 8MM
	Ball Electrode Available Size • Ø 2.3MM • Ø 4.3MM		Needle Electrode Available Size • Ø 4MM • Ø 7MM • Ø 10MM		Blade Electrode Available Size • Ø 25MM

Micro Precision. Maximum Performance




Ultra-Precision Design

- **Equipped with 0.2MM fine tungsten wire tip** for ultra-precise incisions, minimizing scar-ring and thermal spread.
- Ideal for dermatology, cosmetic surgery, dentistry, and ENT applications.



Safety & Sterility

- Gold-plated copper provides excellent **conductivity, corrosion resistance**, and a non-stick surface for optimal performance.
- Pre-sterilized, single-use design ensures safety and minimizes cross-contamination risk.



Flexible & Adaptable

- **Gold-plated copper base** allows for manual bending.
- Adapting easily to anatomical contours for clean, precise cuts with minimal tissue trauma.



Variety of Options

- **Multiple tip options** - loop, ball, needle, and blade - to suit diverse clinical applications.
- Available in $\varnothing 1.6\text{mm}$ and $\varnothing 2.36\text{mm}$ shaft diameters.
- Compatible with most ESU generators for cutting and coagulation.

Microfin Blade Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	WORKING LENGTH(MM)
EM1B2S	Disposable Microfin Blade Electrode	1.6MM	25MM
EM2B2S	Disposable Microfin Blade Electrode	2.36MM	25MM

Microfin Ball Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	BALL DIAMETER(MM)
EM1S2A	Disposable Microfin Ball Electrode	1.6MM	2.3MM
EM2S2A	Disposable Microfin Ball Electrode	2.36MM	4.3MM
EM1S2P	Disposable Microfin Ball Electrode	1.6MM	2.3MM
EM2S2P	Disposable Microfin Ball Electrode	2.36MM	4.3MM

Microfin Bend Needle Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	WORKING LENGTH(MM)
EM1BR2R	Disposable Microfin Bend Needle Electrode	1.6MM	10MM
EM2BR2R	Disposable Microfin Bend Needle Electrode	2.36MM	10MM

Microfin Bend Wire Needle Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	WORKING LENGTH(MM)
EL1BR2R	Disposable Microfin Tungsten Wire Needle Electrode	1.6MM	10MM
EL2BR2R	Disposable Microfin Tungsten Wire Needle Electrode	2.36MM	10MM

Microfin Fine Wire Needle Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	WORKING LENGTH(MM)
EL1WN2R	Disposable Microfin Tungsten Fine Wire Needle Electrode	1.6MM	10MM
EL2WN2R	Disposable Microfin Tungsten Fine Wire Needle Electrode	2.36MM	10MM

Microfin Round Loop Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	LOOP DIAMETER (MM)
EL1R2B	Disposable Microfin Round Loop Tungsten Wire Electrode	1.6MM	4MM
EL2R2B	Disposable Microfin Round Loop Tungsten Wire Electrode	2.36MM	4MM
EL1R2Q	Disposable Microfin Round Loop Tungsten Wire Electrode	1.6MM	7MM
EL2R2Q	Disposable Microfin Round Loop Tungsten Wire Electrode	2.36MM	7MM
EL1R2R	Disposable Microfin Round Loop Tungsten Wire Electrode	1.6MM	10MM
EL2R2R	Disposable Microfin Round Loop Tungsten Wire Electrode	2.36MM	10MM
EL1R2E	Disposable Microfin Round Loop Tungsten Wire Electrode	1.6MM	13MM
EL2R2E	Disposable Microfin Round Loop Tungsten Wire Electrode	2.36MM	13MM

Microfin Diamond Loop Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	LOOP DIAMETER (MM)
EL1D20	Disposable Microfin Diamond Loop Tungsten Wire Electrode	1.6MM	3MM
EL2D20	Disposable Microfin Diamond Loop Tungsten Wire Electrode	2.36MM	3MM
EL1D2C	Disposable Microfin Diamond Loop Tungsten Wire Electrode	1.6MM	5MM
EL2D2C	Disposable Microfin Diamond Loop Tungsten Wire Electrode	2.36MM	5MM
EL1D2D	Disposable Microfin Diamond Loop Tungsten Wire Electrode	1.6MM	8MM
EL2D2D	Disposable Microfin Diamond Loop Tungsten Wire Electrode	2.36MM	8MM

Microfin Triangle Loop Electrode



PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER (MM)	LOOP DIAMETER (MM)
EL1T2B	Disposable Microfin Triangle Loop Tungsten Wire Electrode	1.6MM	4MM
EL2T2B	Disposable Microfin Triangle Loop Tungsten Wire Electrode	2.36MM	4MM
EL1T2Q	Disposable Microfin Triangle Loop Tungsten Wire Electrode	1.6MM	7MM
EL2T2Q	Disposable Microfin Triangle Loop Tungsten Wire Electrode	2.36MM	7MM
EL1T2R	Disposable Microfin Triangle Loop Tungsten Wire Electrode	1.6MM	10MM
EL2T2R	Disposable Microfin Triangle Loop Tungsten Wire Electrode	2.36MM	10MM

Disposable Micro Needle Electrode

Needle Sharp Precision For Fine Electrosurgical Work





Promed's MicroFin Needle Electrode offers a variety of configurations for different surgical procedures. Its ultra-sharp design, combined with fine tungsten wire, ensures precise control and minimizes tissue damage during accurate dissection. The fine needle tip concentrates electrosurgical energy, maintaining sharpness and allowing for lower power settings.



Ø1.6 MM , Ø2.36 MM available

Various Electrode length 50 / 60 / 102 / 152MM

Needle shaft size Ø0.3MM , Ø0.5MM , Ø0.7MM available



Universal Compatibility:

Designed for seamless integration with most standard electrosurgical handpieces.

Ultra-Fine Tungsten Alloy Tip:

Crafted from high heat-resistant tungsten alloy, the Ø0.032MM ultra-sharp, highly polished needle maintains its sharpness even under intense operating conditions.

Enhanced Safety Design:

Equipped with precision shrink wrap insulation to minimize the risk of unintended electrical contact.

Energy Efficient:

The superior conductivity of the tungsten tip allows for effective performance at lower power settings, reducing thermal damage to surrounding tissue

Disposable MicroFin Needle Electrode - Straight Tip

- Tungsten Wire , Ultra Sharp, Straight Tip
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM) LENGTH(MM)			
	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
Post Diameter 2.36MM		PRODUCT CODE	
50MM	EM2305-1	EM2505-1	EM2705-1
60MM	EM2306-1	EM2506-1	EM2706-1
70MM	EM2307-1	EM2507-1	EM2707-1
102MM	EM2310-1	EM2510-1	EM2710-1
152MM	EM2315-1	EM2515-1	EM2715-1
Post Diameter 1.6MM		PRODUCT CODE	
50MM	EM1305-1	EM1505-1	EM1705-1
60MM	EM1306-1	EM1506-1	EM1706-1
70MM	EM1307-1	EM1507-1	EM1707-1
102MM	EM1310-1	EM1510-1	EM1710-1
152MM	EM1315-1	EM1515-1	EM1715-1

Disposable MicroFin Needle Electrode - Straight Tip With Hex Lock

- Tungsten Wire , Ultra Sharp, Straight Tip, Hex Lock
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM) LENGTH(MM)			
	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
Post Diameter 2.36MM		PRODUCT CODE	
102MM	EM2310-1H	EM2510-1H	EM2710-1H
152MM	EM2315-1H	EM2515-1H	EM2715-1H

Disposable MicroFin Needle Electrode - 45° Angled, 3MM

- Tungsten Wire , Ultra Sharp, 45° Angled, 3mm
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM)	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
LENGTH(MM)			
Post Diameter 2.36MM		PRODUCT CODE	
50MM	EM2305-2	EM2505-2	EM2705-2
60MM	EM2306-2	EM2506-2	EM2706-2
70MM	EM2307-2	EM2507-2	EM2707-2
102MM	EM2310-2	EM2510-2	EM2710-2
152MM	EM2315-2	EM2515-2	EM2715-2
Post Diameter 1.6MM		PRODUCT CODE	
50MM	EM1305-2	EM1505-2	EM1705-2
60MM	EM1306-2	EM1506-2	EM1706-2
70MM	EM1307-2	EM1507-2	EM1707-2
102MM	EM1310-2	EM1510-2	EM1710-2
152MM	EM1315-2	EM1515-2	EM1715-2

Disposable MicroFin Needle Electrode - 45° Angled, 3MM With Hex Lock

- Tungsten Wire , Ultra Sharp, 45° Angled, 3mm, Hex Lock
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM)	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
LENGTH(MM)			
Post Diameter 2.36MM		PRODUCT CODE	
102MM	EM2310-2H	EM2510-2H	EM2710-2H
152MM	EM2315-2H	EM2515-2H	EM2715-2H

Disposable MicroFin Needle Electrode - 45° Angled, 10MM

- Tungsten Wire , Ultra Sharp, 45° Angled, 10mm
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM)			
	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
LENGTH(MM)			
Post Diameter 2.36MM		PRODUCT CODE	
60MM	EM2306-3	EM2506-3	EM2706-3
70MM	EM2307-3	EM2507-3	EM2707-3
102MM	EM2310-3	EM2510-3	EM2710-3
152MM	EM2315-3	EM2515-3	EM2715-3
Post Diameter 1.6mm		PRODUCT CODE	
60MM	EM1306-3	EM1506-3	EM1706-3
70MM	EM1307-3	EM1507-3	EM1707-3
102MM	EM1310-3	EM1510-3	EM1710-3
152MM	EM1315-3	EM1515-3	EM1715-3

Disposable MicroFin Needle Electrode - 45° Angled, 10MM With Hex Lock

- Tungsten Wire , Ultra Sharp, 45° Angled, 10mm, Hex Lock
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM)			
	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
LENGTH(MM)			
Post Diameter 2.36MM		PRODUCT CODE	
102MM	EM2310-3H	EM2510-3H	EM2710-3H
152MM	EM2315-3H	EM2515-3H	EM2715-3H

Disposable MicroFin Needle Electrode - 90° Angled, 10MM

- Tungsten Wire , Ultra Sharp, 90° Angled, 10mm
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM)			
	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
LENGTH(MM)			
Post Diameter 2.36MM		PRODUCT CODE	
60MM	EM2306-4	EM2506-4	EM2706-4
70MM	EM2307-4	EM2507-4	EM2707-4
102MM	EM2310-4	EM2510-4	EM2710-4
152MM	EM2315-4	EM2515-4	EM2715-4
Post Diameter 1.6MM		PRODUCT CODE	
60MM	EM1306-4	EM1506-4	EM1706-4
70MM	EM1307-4	EM1507-1	EM1707-4
102MM	EM1310-4	EM1510-1	EM1710-4
152MM	EM1315-4	EM1515-1	EM1715-4

Disposable MicroFin Needle Electrode - 90° Angled, 10MM With Hex Lock

- Tungsten Wire , Ultra Sharp, 90° Angled, 10mm, Hex Lock
- Sterilized Package, 10pcs per Box



NEEDLE SHAFT DIAMETER(MM)			
	Ø 0.3MM	Ø 0.5MM	Ø 0.7MM
LENGTH(MM)			
Post Diameter 2.36MM		PRODUCT CODE	
102MM	EM2310-4H	EM2510-4H	EM2710-4H
152MM	EM2315-4H	EM2515-4H	EM2715-4H

Disposable LLETZ Electrode

Optimized For Cervical Procedure



Promed's gynecology loop electrode ensures precise excision with minimal thermal damage, promoting faster healing and reducing scarring. It delivers consistent energy for accurate cuts, minimizing tissue adhesion and improving procedural efficiency.

DISPOSABLE LLETZ ELECTRODE

Shaft Diameter	2.36MM
Material	Tungsten Wire + Medical-Grade Polymer
Design	Overmold Hex Lock
Shape Available	Loop & Square Shape
Loop Electrode Diameter(MM)	5MM,10MM,15MM,20MM
Square Electrode Diameter(MM)	5MM,8MM,10MM
Packing	1pcs/Pouch,10pcs/Box



FEATURES

Fine Tungsten Wire Construction	High conductivity and durability for precise cutting, minimizing adhesion and ensuring consistent energy.
Precision Design	2.36mm shaft for optimal maneuverability; sterile and single-use to reduce contamination risks.
Overmold Shaft with HEX Lock	Ergonomic, non-slip grip for enhanced control during procedures
Multiple Sizes & Shapes	Available in various loop diameters (e.g., 5mm, 10mm, 15mm, 20mm) for excision, and square electrodes for coagulation.



EL2V64, Loop (A)5MM x (B)5MM

PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER
EL2V64	LLETZ Loop 5*5 MM Electrode	Ø2.36MM



EL2V60, Loop (A)10MM x (B)7MM

PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER
EL2V61	LLETZ Loop 10*5MM Electrode	Ø2.36MM
EL2V60	LLETZ Loop 10*7MM Electrode	Ø2.36MM
EL2V65	LLETZ Loop 10*10MM Electrode	Ø2.36MM



EL2V6B, Loop (A)15MM x (B)5MM

PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER
EL2V6B	LLETZ Loop 15*5MM Electrode	Ø2.36MM
EL2V6C	LLETZ Loop 15*8MM Electrode	Ø2.36MM
EL2V6P	LLETZ Loop 15*10MM Electrode	Ø2.36MM
EL2V6J	LLETZ Loop 15*12MM Electrode	Ø2.36MM
EL2V6K	LLETZ Loop 15*15MM Electrode	Ø2.36MM



EL2V6F, Loop (A)20MM x (B)15MM

PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER
EL2V6D	LLETZ Loop 20*8MM Electrode	Ø2.36MM
EL2V6E	LLETZ Loop 20*10MM Electrode	Ø2.36MM
EL2V6H	LLETZ Loop 20*12MM Electrode	Ø2.36MM
EL2V6F	LLETZ Loop 20*15MM Electrode	Ø2.36MM
EL2V6G	LLETZ Loop 20*20MM Electrode	Ø2.36MM



EL2V74, Square (A)5MM x (B)5MM

PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER
EL2V74	LLETZ Square 5*5MM Electrode	Ø2.36MM



EL2V7A, Loop Square (A)10MM x (B)8MM

PRODUCT CODE	PRODUCT DESCRIPTION	SHAFT DIAMETER
EL2V71	LLETZ Square 10*5MM Electrode	Ø2.36MM
EL2V7A	LLETZ Square 10*8MM Electrode	Ø2.36MM
EL2V75	LLETZ Square 10*10MM Electrode	Ø2.36MM




Disposable Laparoscopic Electrode

Reliable Energy For Minimally Invasive Surgery



Promed Supply offers high-quality disposable monopolar laparoscopic electrodes for enhanced safety in minimally invasive surgeries. With various tip styles and stainless steel or Teflon coatings, they provide versatile options for surgeons and fit standard electrosurgery pencils.



- **Shaft diameter:** $\varnothing$ 2.36mm or $\varnothing$ 4mm, fit most electrosurgery pencils
- **Electrode diameter:** $\varnothing$ 5mm
- **Length:** 33cm or 45cm
- **Tip styles:** L-Wire, L-Hook, J-Wire, J-Hook, Needle, Blade, Ball, etc.
- **Finishing:** Stainless Steel or Non Stick PTFE Teflon Coated



Optional Finishing



Stainless Steel



Teflon Coated

33cm, Ø2.36mm Shaft, L Wire Electrode

45cm, Ø4mm Shaft, Needle Electrode

Disposable Laparoscopic Electrode, 33cm

- Shaft Length 33CM
- Packaged Sterile, 6pcs per Box



PRODUCT	TIP STYLE	SHAFT DIAMETER (MM)	SHAFT LENGTH -33MM	
			OPTIONAL FINISHING	
			STAINLESS	TEFLON COAT
	L- Wire	2.36MM	LE-LW33S2	LE-LW33T2
		4.0MM	LE-LW33S4	LE-LW33T4
	L-Hook	2.36MM	LE-LH33S2	LE-LH33T2
		4.0MM	LE-LH33S4	LE-LH33T4
	J-Wire	2.36MM	LE-JW33S2	LE-JW33T2
		4.0MM	LE-JW33S4	LE-JW33T4
	J-Hook	2.36MM	LE-JH33S2	LE-JH33T2
		4.0MM	LE-JH33S4	LE-JH33T4
	Round Wire	2.36MM	LE-RW33S2	LE-RW33T2
		4.0MM	LE-RW33S4	LE-RW33T4
	Ball	2.36MM	LE-BA33S2	LE-BA33T2
		4.0MM	LE-BA33S4	LE-BA33T4
	Needle	2.36MM	LE-NE33S2	LE-NE33T2
		4.0MM	LE-NE33S4	LE-NE33T4
	Concial Needle Tip	2.36MM	LE-CN33S2	LE-CN33T2
		4.0MM	LE-CN33S4	LE-CN33T4
	Blade	2.36MM	LE-BL33S2	LE-BL33T2
		4.0MM	LE-BL33S4	LE-BL33T4
	Curve Spatula	2.36MM	LE-SP33S2	LE-SP33T2
		4.0MM	LE-SP33S4	LE-SP33T4
	Round Spatula	2.36MM	LE-RS33S2	LE-RS33T2
		4.0MM	LE-RS33S4	LE-RS33T4

Disposable Laparoscopic Electrode, 45cm

- Shaft Length 45CM
- Packaged Sterile, 6pcs per Box



PRODUCT	TIP STYLE	SHAFT DIAMETER (MM)	SHAFT LENGTH -45MM	
			OPTIONAL FINISHING	
			STAINLESS	TEFLON COAT
	L- Wire	ø2.36MM	LE-LW45S2	LE-LW45T2
		ø4.0MM	LE-LW45S4	LE-LW45T4
	L-Hook	ø2.36MM	LE-LH45S2	LE-LH45T2
		ø4.0MM	LE-LH45S4	LE-LH45T4
	J-Wire	ø2.36MM	LE-JW45S2	LE-JW45T2
		ø4.0MM	LE-JW45S4	LE-JW45T4
	J-Hook	ø2.36MM	LE-JH45S2	LE-JH45T2
		ø4.0MM	LE-JH45S4	LE-JH45T4
	Round Wire	ø2.36MM	LE-RW45S2	LE-RW45T2
		ø4.0MM	LE-RW45S4	LE-RW45T4
	Ball	ø2.36MM	LE-BA45S2	LE-BA45T2
		ø4.0MM	LE-BA45S4	LE-BA45T4
	Needle	ø2.36MM	LE-NE45S2	LE-NE45T2
		ø4.0MM	LE-NE45S4	LE-NE45T4
	Concise Needle Tip	ø2.36MM	LE-CN45S2	LE-CN45T2
		ø4.0MM	LE-CN45S4	LE-CN45T4
	Blade	ø2.36MM	LE-BL45S2	LE-BL45T2
		ø4.0MM	LE-BL45S4	LE-BL45T4
	Curve Spatula	ø2.36MM	LE-SP45S2	LE-SP45T2
		ø4.0MM	LE-SP45S4	LE-SP45T4
	Round Spatula	ø2.36MM	LE-RS45S2	LE-RS45T2
		ø4.0MM	LE-RS45S4	LE-RS45T4

Disposable Arthroscopic Electrode

Flexible & Focused energy for Arthroscopy





Promed Arthroscopic Electrodes are designed for cutting and coagulating soft tissue in arthroscopic procedures. Compatible with $\varnothing 2.36\text{MM}$ shaft electro-surgical pencils, they are disposable and available in various tip configurations for use with irrigating solutions.

Insulation extends to the electrode tip, ensuring precise and controlled energy delivery.

The transparent insulation enhances the surgeon's visibility of the electrode tip for improved precision.

The stabilizing cap ensures a secure, safe connection between the electrode and pencil operation.

Compatible with Promed's PD702 pencil for easy conversion.



90° Hook Electrode, Insulated



45° Hook Electrode, Insulated



Angled Blade Electrode, Insulated



Angled Blade Electrode

Optional Anti-Rotation




with hexagon lock




without hexagon lock

Hexagon lock to pair with
PROMED PD702 ES Pencil

Optional Finishing



Stainless Steel



Teflon Coated

Disposable Arthroscopic Electrode

- Shaft Diameter: $\varnothing$ 2.36MM; Length: 175MM (cap to electrode)
- Available in the finishing of Stainless Steel or Teflon Coating
- Sterilization Packaging
 - Electrodes Only :1pcs/Pouch, 50pcs/Box
 - PD702 Pencil: 1pcs/Tyvek Pouch, 10pcs/Box (Optional)



PRODUCT	PRODUCT DESCRIPTION	PRODUCT CODE			
		STAINLESS STEEL TIP	TEFLON COAT TIP	STAINLESS STEEL TIP	TEFLON COAT TIP
Options -W or W/ Hexagon Lock		 without hexagon lock		 with hexagon lock	
	Transparent Insulated Hook Electrode, 90°	EA801S	EA801T	EA801S-H	EA801T-H
	Transparent Insulated Hook Electrode, 45°	EA802S	EA802T	EA802S-H	EA802T-H
	Transparent Insulated Angled Blade Electrode	EA803S	EA803T	EA803S-H	EA803T-H
	Uninsulated Angled Blade Electrode	EA804S	EA804T	EA804S-H	EA804T-H

PRODUCT		PRODUCT DESCRIPTION	STAINLESS STEEL TIP	TEFLON COAT TIP
		PS702 Pencil with Transparent Insulated Hook Electrode, 90°	PDA801S-H	PDA801T-H
		PS702 Pencil with Transparent Insulated Hook Electrode, 45°	PDA802S-H	PDA802T-H
		PS702 Pencil with Transparent Insulated Angled Blade Electrode	PDA803S-H	PDA803T-H
		PS702 Pencil with Uninsulated Angled Blade Electrode	PDA804S-H	PDA804T-H