


Disposable Laparoscopic Electrode

Reliable Energy For Minimally Invasive Surgery



Promed Supply offers high-quality disposable monopolar laparoscopic electrodes for enhanced safety in minimally invasive surgeries. With various tip styles and stainless steel or Teflon coatings, they provide versatile options for surgeons and fit standard electrosurgery pencils.



- **Shaft diameter:** \varnothing 2.36mm or \varnothing 4mm, fit most electrosurgery pencils
- **Electrode diameter:** \varnothing 5mm
- **Length:** 33cm or 45cm
- **Tip styles:** L-Wire, L-Hook, J-Wire, J-Hook, Needle, Blade, Ball, etc.
- **Finishing:** Stainless Steel or Non Stick PTFE Teflon Coated



Optional Finishing



Stainless Steel



Teflon Coated












33cm, Ø2.36mm Shaft, L Wire Electrode

45cm, Ø4mm Shaft, Needle Electrode

Disposable Laparoscopic Electrode, 33cm

- Shaft Length 33CM
- Packaged Sterile, 6pcs per Box














PRODUCT	TIP STYLE	SHAFT DIAMETER (MM)	SHAFT LENGTH -33MM	
			OPTIONAL FINISHING	
			STAINLESS	TEFLON COAT
	L- Wire	2.36MM	LE-LW33S2	LE-LW33T2
		4.0MM	LE-LW33S4	LE-LW33T4
	L-Hook	2.36MM	LE-LH33S2	LE-LH33T2
		4.0MM	LE-LH33S4	LE-LH33T4
	J-Wire	2.36MM	LE-JW33S2	LE-JW33T2
		4.0MM	LE-JW33S4	LE-JW33T4
	J-Hook	2.36MM	LE-JH33S2	LE-JH33T2
		4.0MM	LE-JH33S4	LE-JH33T4
	Round Wire	2.36MM	LE-RW33S2	LE-RW33T2
		4.0MM	LE-RW33S4	LE-RW33T4
	Ball	2.36MM	LE-BA33S2	LE-BA33T2
		4.0MM	LE-BA33S4	LE-BA33T4
	Needle	2.36MM	LE-NE33S2	LE-NE33T2
		4.0MM	LE-NE33S4	LE-NE33T4
	Concial Needle Tip	2.36MM	LE-CN33S2	LE-CN33T2
		4.0MM	LE-CN33S4	LE-CN33T4
	Blade	2.36MM	LE-BL33S2	LE-BL33T2
		4.0MM	LE-BL33S4	LE-BL33T4
	Curve Spatula	2.36MM	LE-SP33S2	LE-SP33T2
		4.0MM	LE-SP33S4	LE-SP33T4
	Round Spatula	2.36MM	LE-RS33S2	LE-RS33T2
		4.0MM	LE-RS33S4	LE-RS33T4

Disposable Laparoscopic Electrode, 45cm

- Shaft Length 45CM
- Packaged Sterile, 6pcs per Box



PRODUCT	TIP STYLE	SHAFT DIAMETER (MM)	SHAFT LENGTH -45MM	
			OPTIONAL FINISHING	
			STAINLESS	TEFLON COAT
	L- Wire	2.36MM	LE-LW45S2	LE-LW45T2
		4.0MM	LE-LW45S4	LE-LW45T4
	L-Hook	2.36MM	LE-LH45S2	LE-LH45T2
		4.0MM	LE-LH45S4	LE-LH45T4
	J-Wire	2.36MM	LE-JW45S2	LE-JW45T2
		4.0MM	LE-JW45S4	LE-JW45T4
	J-Hook	2.36MM	LE-JH45S2	LE-JH45T2
		4.0MM	LE-JH45S4	LE-JH45T4
	Round Wire	2.36MM	LE-RW45S2	LE-RW45T2
		4.0MM	LE-RW45S4	LE-RW45T4
	Ball	2.36MM	LE-BA45S2	LE-BA45T2
		4.0MM	LE-BA45S4	LE-BA45T4
	Needle	2.36MM	LE-NE45S2	LE-NE45T2
		4.0MM	LE-NE45S4	LE-NE45T4
	Concise Needle Tip	2.36MM	LE-CN45S2	LE-CN45T2
		4.0MM	LE-CN45S4	LE-CN45T4
	Blade	2.36MM	LE-BL45S2	LE-BL45T2
		4.0MM	LE-BL45S4	LE-BL45T4
	Curve Spatula	2.36MM	LE-SP45S2	LE-SP45T2
		4.0MM	LE-SP45S4	LE-SP45T4
	Round Spatula	2.36MM	LE-RS45S2	LE-RS45T2
		4.0MM	LE-RS45S4	LE-RS45T4