



Disposable Laparoscopic Instrument

High-Performance Tool For MIS



Promed's disposable monopolar laparoscopic instruments are specifically engineered for minimally invasive surgical procedures, prioritizing surgeon comfort, safety, and ease of use. Available with either single or double action jaws, they deliver precise control, enhanced tactile feedback, and are equipped with a standard monopolar connection for seamless integration with electrosurgical systems.
















- 1 Ergonomic handles for optimal comfort and control
- 2 Large finger rings to reduce fatigue during prolonged use
- 3 360° fingertip rotation for precise tip orientation
- 4 On/Off ratchet mechanism for enhanced precision and control



- 5 Smooth jaw movement with maximum opening
- 6 Sterilized and ready-to-use for OR staff
- 7 330mm shaft length
- 8 Sterile Packing
- 9 10pcs/Box



PRODUCT	STYLE	ARTICULATION	JAW LENGTH (MM)	SHAFT LENGTH (MM)	 WITH RATCHET	 WITHOUT RATCHET
	Straight Scissor	Double	13MM	330MM	/	L09S-330-2
	Curve Scissor	Double	17MM	330MM	/	L09C-330-2
	Hook Scissor	Double	15MM	330MM	/	L09G-330-1
	Maryland Forceps	Double	18MM	330MM	L18ML-1-330-2	L18ML-0-330-2
	Duck Bill Forceps	Double	18MM	330MM	L18DB-1-330-2	L18DB-0-330-2
	Endoclinch Forceps	Double	20MM	330MM	L19EC-1-330-2	L19EC-0-330-2
	Claw Forceps	Single	22MM	330MM	L21-330-1	/
	Mixer Dissector	Double	21MM	330MM	L18M-0-330-2	/
	Badcock Forceps	Double	28MM	330MM	L19BC-1-330-2	/
	Biopsy Forceps	Single	17MM	330MM	/	L21-330-1
	Biopsy Forceps with Pin	Single	17MM	330MM	/	L05-330-2