



# Disposable Laparoscopic Instrument

High-Performance Tool For MIS



Promed's disposable monopolar laparoscopic instruments are specifically engineered for minimally invasive surgical procedures, prioritizing surgeon comfort, safety, and ease of use. Available with either single or double action jaws, they deliver precise control, enhanced tactile feedback, and are equipped with a standard monopolar connection for seamless integration with electrosurgical systems.
















- 1 Ergonomic handles for optimal comfort and control
- 2 Large finger rings to reduce fatigue during prolonged use
- 3 360° fingertip rotation for precise tip orientation
- 4 On/Off ratchet mechanism for enhanced precision and control



- 5 Smooth jaw movement with maximum opening
- 6 Sterilized and ready-to-use for OR staff
- 7 330mm shaft length
- 8 Sterile Packing
- 9 10pcs/Box




PRODUCT	STYLE	ARTICULATION	JAW LENGTH (MM)	SHAFT LENGTH (MM)	 WITH RATCHET	 WITHOUT RATCHET
	<b>Straight Scissor</b>	Double	13MM	330MM	/	L09S-330-2
	<b>Curve Scissor</b>	Double	17MM	330MM	/	L09C-330-2
	<b>Hook Scissor</b>	Double	15MM	330MM	/	L09G-330-1
	<b>Maryland Forceps</b>	Double	18MM	330MM	L18ML-1-330-2	L18ML-0-330-2
	<b>Duck Bill Forceps</b>	Double	18MM	330MM	L18DB-1-330-2	L18DB-0-330-2
	<b>Endoclinch Forceps</b>	Double	20MM	330MM	L19EC-1-330-2	L19EC-0-330-2
	<b>Claw Forceps</b>	Single	22MM	330MM	L21-330-1	/
	<b>Mixer Dissector</b>	Double	21MM	330MM	L18M-0-330-2	/
	<b>Badcock Forceps</b>	Double	28MM	330MM	L19BC-1-330-2	/
	<b>Biopsy Forceps</b>	Single	17MM	330MM	/	L21-330-1
	<b>Biopsy Forceps with Pin</b>	Single	17MM	330MM	/	L05-330-2

# Suction Coagulator

Dual Function :Suction & Coagulation



Tube Diameter 8Fr/10Fr/12Fr



Type	Footswitch or Hand Switch
Size Available	8Fr,10Fr,12Fr
Suction:	Place Finger Over Suction Vent on Handle
Handle:	Ergonomically Designed White ABS
Cannula Insulation:	Blue Polyolefin
Obturator:	Polypropylene Cleaning Tube and Protective Insert for Bending
Suction Tubing:	6mm or 5mm Tubing Compatible
Cable:	3m(10ft)/5m(15ft)for both Footswitch and Handswitch

## Handswitch Suction Coagulators

PRODUCT	CODE	TUBE DIAMETER	CABLE LENGTH(M)	PACKING
	SC08H	8fr	3M	15/Case
	SC10H	10fr	3M	15/Case
	SC12H	12fr	3M	15/Case

## Footswitch Suction Coagulators

PRODUCT	CODE	TUBE DIAMETER	CABLE LENGTH(M)	PACKING
	SC08F	8fr	3M	15/Case
	SC10F	10fr	3M	15/Case
	SC12F	12fr	3M	15/Case




# Disposable Suction Irrigator

*Clean Precision , Seamless Control*

Engineered to streamline surgical workflows, Promed's Suction Irrigator integrates simultaneous suction and irrigation in a unified system. Featuring an ergonomic design and standardized connectivity, it maintains optimal visibility during procedures while supporting precision in fluid evacuation and targeted irrigation.

GENERAL SURGERY SOLUTION



## Order Information

PRODUCT CODE	PRODUCT DESCRIPTION	PACKING
SI0528PD	Disposable Suction Irrigator, $\varnothing$ 5.36mmx28cm Composites Probe, Double Spikes	20/Case
SI0533PD	Disposable Suction Irrigator, $\varnothing$ 5.36mmx 33cm Composites Probe, Double Spikes	20/Case
S10545PD	Disposable Suction Irrigator, $\varnothing$ 5.36mmx45cm Composites Probe, Double Spikes	20/Case
S10528PS	Disposable Suction Irrigator, $\varnothing$ 5.36mm X28cm Composites Probe, Single Spike	20/Case
S10533PS	Disposable Suction Irrigator, $\varnothing$ 5.36mmx 33cm Composites Probe, Single Spike	20/Case
SI0545PS	Disposable Suction Irrigator, $\varnothing$ 5.36mm X45cm Composites Probe, Single Spike	20/Case