

# Disposable Microfin Electrosurgical Pencil

Exceptional Precision for Micro-Electrosurgery



Cut and coagulation buttons are conveniently positioned near the surgical site for improved control, precision, and smooth electrode handling.

Compatible with ultra-fine  $\varnothing 1.6\text{mm}$  electrode for superior surgical precision.

Lightweight, compact design for easy handling and reduced fatigue.

**PD833W -**  
 $\varnothing 1.6\text{MM}$

**PD833 -**  
 $\varnothing 2.36\text{MM}$

Compatible with standard  $\varnothing 2.36\text{mm}$  electrodes for versatile surgical applications.

\*\*\* Electrode is not included with the ESU pencil \*\*\*

## Disposable Microfin Electrosurgical Pencil

- Sterile Packing
- Packing: 10 pcs/Box



PRODUCT CODE	ELECTRODE DIAMETER	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
PD833W	$\varnothing 1.6\text{MM}$	106MM	3M	3 Prong Plug
PD833	$\varnothing 2.36\text{MM}$	106MM	3M	3 Prong Plug