

Disposable Skin Stapler & Remover

Fast & Safe & Consistent



Disposable Skin Stapler

- Precise staple forming ensures optimal wound approximation and healing.
- Consistent, smooth, and safe actuation provides accurate staple depth.
- Preloaded with medical-grade stainless steel staples for reliable closure.
- EO is sterilized and designed for single-use to guarantee maximum hygiene.
- A transparent window allows easy monitoring of remaining staples.

Staple Dimensions

Open	14mm  3.2mm
Closed	W 6.5+0.5/-0.2mm  H 4.7+0.5/-0.4mm



PRODUCT CODE	PRODUCT DESCRIPTION	Staple	STERILE PACKING
SW904-15W	Disposable Skin Stapler SW904,15W	15PCS	6/Case
SW904-25W	Disposable Skin Stapler SW904,25W	25PCS	6/Case
SW904-35W	Disposable Skin Stapler SW904,35W	35PCS	6/Case

Disposable Skin Staple Remover

Promed's Disposable Skin Stapler Remover features an angled tip for easy access and precise staple removal. Designed for gentle tissue handling, it minimizes patient discomfort during procedures. Its simple, user-friendly mechanism ensures efficient and safe use by medical professionals. Sterile and disposable, it supports optimal hygiene and convenience in clinical settings.



Precision-Engineered Tip	Designed to securely grip and remove staples without causing trauma to surrounding tissue.
Angled Tip Design	Allows for easy access and optimal visibility, especially in hard-to-reach areas.
Ergonomic Handle	Comfortable, non-slip grip for controlled and accurate removal.
Sterile & Single-Use	Pre-sterilized and disposable to prevent cross-contamination and ensure patient safety.
Lightweight & Easy to Use	Compact and lightweight design minimizes hand fatigue and enhances ease of use during procedures.
Durable Medical-Grade Material	Made from high-quality plastic and stainless steel components for reliable performance.

PRODUCT CODE	PRODUCT DESCRIPTION	Staple	STERILE PACKING
SR00D1	Disposable Skin Staple Remover	15PCS	10/Case