



# Reusable Electrosurgical Pencil

Reliable, Reusable, and Cost-Effective



GENERAL SURGERY SOLUTION

Reusable for up to 30/50/100 cycles, reducing long-term operating costs.

Push-button activation enables precise and efficient electrosurgical performance.

Ergonomic and lightweight design for enhanced surgeon comfort and control.

Compatible with 2.36mm electrodes and most standard electrosurgical units.

## Reusable Electrosurgical Pencil

- Fit with  $\varnothing$  2.36mm Electrode
- Sterile Packing
- Packing: 10pcs/Box

134°C  
Autoclave



PRODUCT CODE	BLADE TYPE	REUSABLE	HOUSING LENGTH(MM)	CABLE LENGTH(M)	PLUG TYPE
PR306P	Standard Reusable Blade	30 Cycles	163MM	3M	3 Prong Plug
PR307P	Standard Reusable Blade	50 Cycles	163MM	3M	3 Prong Plug
PR310P	Standard Reusable Blade	100 Cycles	163MM	3M	3 Prong Plug