

Reusable Bipolar Forceps

Durable Precision For Bipolar Coagulation



Promed's Reusable Bipolar Forceps are available in a variety of lengths, tip sizes, and designs, including Bayonet, Jewelers, Adson, and Straight, to suit diverse surgical needs.

Engineered for optimal performance, the polished stainless-steel tips reduce eschar buildup, ensuring smooth operation. The ergonomic, balanced handles minimize hand fatigue, enhancing comfort and precision. Designed for durability, they withstand up to 20 auto-clave cycles.



Polished Stainless Steel Tips

Secure grasping and coagulation at lower power, reducing tissue damage and eschar buildup.



Insulated Stainless Steel Body

Prevents electrical leakage for enhanced safety.



Precise Handle Spring Tension

Ensuring consistent dissection pressure and maintaining accurate tip alignment.




Platform-Like Tip Design

Increases resistance and enhances coagulation with controlled heat.



Ergonomic, Well-Balanced Handles


Provides balance and comfort reducing hand fatigue.



STYLE	LENGTH(MM)	TIP WIDTH (MM)
Bayonet Straight	180/200/220/240MM	0.5/1.0/2.0MM
Bayonet Angled Down	200/220/240MM	0.5/1.0/2.0MM
Bayonet Angled Up	200/220/240MM	0.5/1.0/2.0MM
Straight	140/180/200/220/240MM	0.5/1.0/2.0MM
Adson	130MM	0.5/1.0MM
Jeweler	110MM	0.5MM

Straight Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack



FST0205, Length 200mm, Straight, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
Straight Tip, 2 Pins Connector					
0.5MM	FST0145	FST0185	FST0205	/	/
1.0MM	FST0141	FST0181	FST0201	FST0221	FST0241
2.0MM	FST0142	FST0182	FST0202	FST0222	FST0242
Straight Tip, Flat Connector					
0.5MM	FST1145	FST1185	FST1205	/	/
1.0MM	FST1141	FST1181	FST1201	FST1221	FST1241
2.0MM	FST1142	FST1182	FST1202	FST1222	FST1242




FAG0205, Length 200mm, Straight, Angled Tip, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
Straight, Angled Tip, 2 Pins Connector					
0.5MM	/	FAG0185	FAG0205	/	/
1.0MM	/	FAG0181	FAG0201	FAG0221	FAG0241
2.0MM	/	FAG0182	FAG0202	FAG0222	FAG0242
Straight, Angled Tip, Flat Connector					
0.5MM	/	FAG1185	FAG1205	/	/
1.0MM	/	FAG1181	FAG1201	FAG1221	FAG1241
2.0MM	/	FAG1182	FAG1202	FAG1222	FAG1242

Bayonet Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack



FBD0205, Length 200mm, Bayonet, Angled Down, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
Bayonet, Angled Down, 2 Pins Connector					
0.5MM	/	/	FBD0205	FBD0225	/
1.0MM	/	/	FBD0201	FBD0221	FBD0241
2.0MM	/	/	FBD0202	FBD0222	/
Bayonet, Angled Down, Flat Connector					
0.5MM	/	/	FBD1205	FBD1225	/
1.0MM	/	/	FBD1201	FBD1221	FBD1241
2.0MM	/	/	FBD1202	FBD1222	/



FBS0205, Length 200mm, Bayonet, Straight, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
Bayonet, Straight Tip, 2 Pins Connector					
0.5MM	/	FBS0185	FBS0205	FBS0225	FBS0245
1.0MM	/	FBS0181	FBS0201	FBS0221	FBS0241
2.0MM	/	FBS0182	FBS0202	FBS0222	/
Bayonet, Straight Tip, Flat Connector					
0.5MM	/	FBS1185	FBS1205	FBS1225	FBS1245
1.0MM	/	FBS1181	FBS1201	FBS1221	FBS1241
2.0MM	/	FBS1182	FBS1202	FBS1222	/




FBU0205, Length 200mm, Bayonet, Angled Up, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)				
	140MM	180MM	200MM	220MM	240MM
Bayonet, Angled Up, 2 Pins Connector					
0.5MM	/	/	FBU0205	FBU0225	/
1.0MM	/	/	FBU0201	FBU0221	FBU0241
2.0MM	/	/	FBU0202	FBU0222	/
Bayonet, Angled Up, Flat Connector					
0.5MM	/	/	FBU1205	FBU1225	/
1.0MM	/	/	FBU1201	FBU1221	FBU1241
2.0MM	/	/	FBU1202	FBU1222	/

Jewelers

Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack




MJE0115, Length 110mm, Jewelers, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)
	110MM
Jeweler, 2 Pins Connector	
0.5MM	MJE0115
1.0MM	/
Jeweler, Flat Connector	
0.5MM	MJE1115
1.0MM	/

Adson

Reusable Bipolar Forceps

- High polished stainless-steel surface
- Packaged non- sterile, 1 each per pack



MAD0135, Length 130mm, Adson, 0.5mm Tip Width, 2 Pins Connector

TIP WIDTH (MM)	FORCEPS LENGTH(MM)
	130MM
Adson, 2 Pins Connector	
0.5MM	MAD0135
1.0MM	MAD0131
Adson, Flat Connector	
0.5MM	MAD1135
1.0MM	MAD1131